NSM ST NICAL
SMART MICRO ST

NSM ST NICAL SMART MICRO ST

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Power supply and PC connector cables are located on the back.

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LED indicator displaying steps of operation.

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Slot for inserting image plate. Eliminates contact to light and fingers.

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Light cover

Hygienic bag

Nical continues in the growth of its own range of phosphor plate reader presenting the innovative model ST with an enlarged field of application and an even better image quality.

The simple and minimalist design of the appliance suits well a modern dental practice thanks to the compact overall dimensions and constitutes the last and most modern project in the sector thanks to a compact and rational electronic 100% SMD, useful to pursuit a drastic reduction of failures and machine stoppage.

The phosphor plate is inserted in the upper part of the machine and the process of reading, erasing and ejection as well as the plate dimension detection is totally automatic.

Any malfunction is monitored by a color code reproduced from a smart led display.

The plate is managed through the light cover, which provides protection from light and ensures cleanliness by avoiding contact with the operator's hands.

The maximization of user-friendliness advocated by NICAL is found again in the imaging software that does not require special adjustment to be made on the image capturing, thanks both to the particular preprocessing, previously carried out at the hardware level during the reading process, and to the automatic adjustment of the exposure, which grants to eliminate over and under exposure. Beside that, the NEW imaging software, based on a multi standard platform which can mix different operating systems and different hardware (either windows or android), offers a very high level of flexibility, system integration and data exchanging never seen before.

The powerful Imaging software offers DICOM interface with RIS and PACS; furthermore, it can be also linked to third-parties practice management application so as to share patient's database and multi printing procedures are allowed.

All of this, together with a superb spatial and wide contrast resolution, characteristic of the unique

technology of the Smart Micro line, makes the iconography of the ST a reference for the market, not only in the middle price segment to which the product belongs to but also to the most qualify high level product segment.

The concept of quality and user-friendliness must get along with that of easy and efficient maintenance during the operative life of the unit.

The Smart Micro ST offers a wide range of image formats such as:

size0 31x22

size1 40x24

size2 41x31

size3 54x27

The four image sizes are automatically detected and identified by the system during the reading process, maintaining the same resolution per mm from size 0 to size size 5.

The plate of the Smart Micro ST are patented by NICAL and they are assembled on the holding, coupling a high sensibility phosphor plate with a flexible magnetic support. This technology offers a longer operative life, an extra radiation shield for the patient and a clever scattering reduction for a better image quality.

Furthermore, the magnetic support allows to move the plate inside the machine during the reading procedure without any mechanic clutch or rubbing, which grants to preserve the phosphor surface integrity for a longer time.

The connection of the unit to the PC or into a network is granted by a double modality, either USB or Ethernet. Two modules, easy to install, are available to enable the USB or the Ethernet connection.

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NICAL SMART MICRO ST technical and performance data

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Caratteristiche PC e Software PC and Software features

Processore PC <i>Processor PC</i>	32/64 bit Intel Dual core
RAM <i>RAM</i>	1GB
Matrice display <i>Matrice display</i>	1280x1024 minimum
Hard Disk <i>Hard Disk</i>	>500 GB
Sistema operativo <i>Windows® operating system</i>	Windows: 10, 8.1, 7, Vista, Mac iOS, Android, opz. per interfaccia Ethernet

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
The light cover and the disposable hygienic bag are available for all sizes used as follows.

<p>SIZE 0 31x22 mm 726x1024 pixel</p>	<p>SIZE 1 40x24 mm 792x1321 pixel</p>	<p>SIZE 3 54x27 mm 891x1783 pixel</p>	<p>SIZE 2 41x31 mm 1024x1354 pixel</p>
---	---	---	--

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The system comes with all the accessories and a complete image plate starter kit allowing immediate use of the NICAL SMART MICRO.

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
Different kinds of acquisition mode

- **FMX**
- **Introrale**
- **OPG***
- **3D***
- **Ceph***
- **Videocamera**
- **New generation TWAIN**
- **Photocamera**
- **Scanner**
- **Import images**

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