



- Clinical Imagery







- Sharp, high resolution images with a 0.4 focal spot

Source focal spot Source focal spot





Eigilieetii

Effect of focal spot size



## Designed with safety in mind



Built-in leaded double shield to protect the operator from leakage radiation



Locking system to prevent unintentional exposure

## **Specifications**

Model	HyperLight	Weight	1.8 kg
Focal Spot	0.4	Battery	2500 mAh x 4
Tube Voltage	65 kV	Exposure Time Range	0.02-2 Sec
Tube Current	2.5 mA	Line voltage	100V-240V, 50Hz/60Hz

HyperLight


44







## Simplify your digital workflow


- Compact size and lightweight can make you operate it by one hand only
- Intuitive user interface











## Maximize your diagnostic potential

Perfect combination



HyperLight

46