

## APEX LOCATOR DPEX V

- ✦ 6<sup>th</sup> generation apex locator for measurement in presence of blood, pus, residual pulp or saliva etc.
- ✦ 97.71% accurate than the conventional radiography method
- ✦ Multi-frequency impedance technology precisely measure irrespective of the tooth type
- ✦ File clip, hook, file probe can be easily Autoclavable to avoid cross-contamination
- ✦ Display: Multi-color 4.5" LCD display
- ✦ Main unit weight: 385gm
- ✦ Dimensions: 138mm x 77mm x 107mm

