


DB686 NanoLED curing light

Small, just right;
Use, just handy.

Small body, powerful function
Light curing & caries detection

2 in 1



Model: DB686 NANO

LED Curing Lights

Order code: DB686 Nano

Technical:

Adapter: AC100-240V, DC5V/1.5A, 50/60Hz

Light curing mode:

Wavelength: 380-520nm

Curing time: 10s/20s


Curable depth: 8mm

Caries detection mode:


Wavelength: 380-420nm



LED light
measure system




Ultra-small all-metal body
outer diameter of only 13.5mm



8 mm deep curing depth
Curing effect into every corner




10mm diameter
large area light emission




Wavelength: 380-520nm
Suitable for most of the resins
in the market

Detection mode:
Simple, Effective
Innovative caries detection



Green: Healthy



Red: Bacterial dentin

Nano can maintain strength even if
it is not focused on the center.
Curing at any angle can also ensure the effect

